

Suture for Orthopedic surgery

FiberPLUS®

Ultra High
MOLECULAR WEIGHT POLYETHYLENE
(UHMWPE)

● FEATURES

- **Strong non-absorbable** suture made from **Ultra High Molecular Weight Polyethylene**
- **Outstanding** straight tensile and knot **strength**
- Flexible, durable
- For **orthopedic** procedures

● SPECIFICATION


- Available size range : USP 2 - 2/0
- Available colors : White, White/Green, White/Blue, White/Violet, White/Black
---> for easier visualization, identification and motion detection.

Needles
AISI
302


Fiber threads

← TENSILE STRENGTH →

Linear Tensile Strength



EP knot Strength



Suture for Orthopedic surgery

FiberPLUS®

- Strong non-absorbable suture made from Ultra High Molecular Weight Polyethylene
- Outstanding straight tensile and knot strength
- Flexible, durable
- For orthopedic procedures

Ultra High
MOLECULAR WEIGHT POLYETHYLENE
(UHMWPE)



**POLYETHYLENE
(UHMWPE)**

(USP 2 - 35")

Color: **White**

(USP 2 - 35")



Color: **White - Green**

(USP 2/0 - 35")

Color: **White**

(USP 2/0 - 35")

Color: **White - Green**

USP Size	Needle's shape and type	Needle size (*)	Length
2	 1/2 circle taper point Ⓞ	25 mm, 30 mm, 35 mm	100 cm
USP Size	Needle's shape and type	Needle size (*)	Length
2/0	 1/2 circle taper point Ⓞ	25 mm, 30 mm, 35 mm	100 cm